RTA®

Plasmid DNA Isolation Kit

Handbook

Date of Issue -- 2011-12

For isolation and purification of plasmid DNA from bacterial samples
For professional use only

09007050 - 50 tests
09007100 - 100 tests

Table of Contents

| | |
|-------------------------------------|-------|
| Kit Contents | 4 |
| Storage and Stability | 4 |
| Intended Use | 5 |
| Product Use Limitations | 5 |
| Introduction | 6-7 |
| Warnings and Precautions | 8 |
| Additional Materials Required | 9 |
| Notes Before Starting | 9 |
| Protocol | 10-11 |
| Troubleshooting Guide | 12 |

Kit Contents

| Materials provided | Quantity | | |
|---------------------------|----------|----------|-----------|
| | 10 tests | 50 tests | 100 tests |
| Solution R | 2 ml | 5.5 ml | 11 ml |
| Solution L | 2 ml | 5.5 ml | 11 ml |
| Solution N | 2 ml | 5.5 ml | 11 ml |
| Solution P | 3 ml | 11 ml | 22 ml |
| Binding Solution | 5 ml | 22 ml | 44 ml |
| Wash Solution* | 2 ml | 8 ml | 16 ml |
| Spin Columns | 10 | 50 | 100 |
| Collection Tubes (2 ml) | 10 | 50 | 100 |
| Collection Tubes (1.5 ml) | 10 | 50 | 100 |
| Handbook | 1 | 1 | 1 |

* For preparation of Wash Solution see Notes Before Starting section.

Storage & Stability

All components of the RTA Plasmid DNA isolation kit except Solution R should be stored at room temperature (15–25°C). It is recommended to store Solution R at 2–8 °C after arrival of the kit for long-term usage. Under these conditions, the kits can be stored for at least 12 months.

Intended Use

RTA Plasmid DNA Isolation Kit is a plasmid extraction system from bacterial cultures. Plasmid DNA is purified from bacterial samples, which can be further used in sensitive molecular as well as diagnostic downstream applications such as PCR, endonuclease digestion, subcloning, transformation, ligation and sequencing reactions.

Product Use Limitations

RTA Plasmid DNA Isolation Kit is intended for general nucleic acid purification purposes and not for any specific organism or for any specific clinical purpose. It is up to the user's discretion whether or not it is suitable for their specific experimental needs.

Introduction


RTA Plasmid DNA Isolation Kit provides a rapid and easy method for purifying highly pure plasmid DNA from bacterial cultures without the use of phenol, chloroform, CsCl, or alcohol precipitation. Nucleic acids bind specifically to the surface of glass or silica materials in the presence of a chaotropic salt. The binding reaction occurs due to the disruption of the organized structure of water molecules and the interaction with the nucleic acids. Since the binding process is specific for nucleic acids, the bound material can be separated and purified from impurities e.g. salts and proteins, by a simple washing step. Nucleic acids elute from the matrix in a low salt buffer or water. The method is divided into three basic steps: preparation of a cleared lysate, adsorption of DNA onto glass microfiber, and washing and subsequent elution of plasmid DNA. The result is purified plasmid DNA that is ready for restriction endonuclease digestion, PCR, subcloning, transformation, ligation, and sequencing reactions.

Plasmid DNA yield can be affected by a number of factors, including plasmid copy number, insert size, the volume of culture processed, culture media and the overall binding capacity of the purification system. Of these factors, plasmid copy number, culture volume and binding capacity have the greatest influence on plasmid yield. Purified plasmid DNA is free of enzymatic inhibitors as measured by restriction endonuclease digestion. The plasmid DNA is predominantly supercoiled with no RNA contamination.

Introduction

(continued)

| | |
|----------|--|
| Yield | Up to 20 µg of plasmid DNA |
| Speed | Takes only 20 minutes to extract plasmid DNA |
| Specific | No detectable genomic DNA or RNA contamination |
| Safe | No phenol/chloroform extraction required |
| Purity | Highly purified plasmid DNA with our specially designed plasmid DNA specific columns |

Warnings and Precautions

Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example; safety glasses, gloves or protective clothing). Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example; fume hood). For additional safety guidelines, consult the MSDS.

The sample-preparation waste can form highly reactive compounds when combined with bleach. DO NOT add bleach or acidic solutions directly to the sample-preparation waste.

Dispose of unused reagents, waste and specimens in accordance with country or local regulations.

Do not pipette by mouth.

Do not eat, drink or smoke in laboratory work areas.

Avoid contact of reagents with the skin, eyes or mucous membranes. If contact does occur, immediately wash with large amounts of water,

Use all pipetting devices and instruments with care and follow the manufacturer's instructions for calibration and quality control; to prevent sample contamination, use new, sterile aerosol barrier or positive displacement DNase-free pipette tips and sterile pipettes.

Handle all materials containing specimens or controls according to Good Laboratory Practices in order to prevent cross-contamination of specimens or controls. Store the kit away from any source of contaminating DNA or RNA, especially amplified nucleic acid.

Do not mix reagents with different lot numbers or substitute reagents from other manufacturers.

Do not use a kit after its expiration date.

Additional Materials Required

- 96-100% ethanol (Molecular biology grade)
- sterilized dH₂O or a low-ionic strength buffer (See the Protocol Step 11)
- vortex
- Microcentrifuge
- Microcentrifuge tubes (1.5 ml)

Notes Before Starting

- Under cool ambient conditions, a precipitate may form in Solution L. If a precipitate is observed in Solution L, incubate the bottle at 37°C for some minutes to dissolve. Always keep bottles tightly closed and the buffers should be at room temperature when performing the protocol.
- Wash Solution is a concentrate. Before using for the first time, add the appropriate amount of ethanol (96-100%) as indicated on the bottle and in the table below:

| <i>Kit Size</i> | <i>Wash Solution (ml)</i> | <i>Ethanol to be added (ml)</i> | <i>Final Volume (ml)</i> |
|-----------------|---------------------------|---------------------------------|--------------------------|
| 10 tests | 2 | 8 | 10 |
| 50 tests | 8 | 32 | 40 |
| 100 tests | 16 | 64 | 80 |

Protocol

- 1. Harvest 1.5 ml of the overnight bacterial culture for 1 minute at 14,000 rpm.**

Discard supernatant.

*Remove all traces of supernatant by inverting the centrifuge tube until all media has been drained. The pellet can then be frozen at -20°C for later use. The protocol is designed for use with 1.5 ml of *E.coli* cell culture. The protocol is flexible and can be scaled to accommodate other sample volumes.*

- 2. Resuspend each sample with 100 μl of Solution R by vigorously vortexing.**

Make sure that the bacterial pellet is completely resuspended and no clumps remain.

- 3. Add 100 μl of Solution L and mix gently by inverting 5 times. Incubate at room temp. for 1 minute.**

DO NOT VORTEX. Vortexing will result in shearing of the genomic DNA, leaving free chromosomal fragments to contaminate the plasmid DNA.

- 4. Add 100 μl of Solution N and mix gently by inverting 5 times. Incubate at room temp. for 3 minutes.**

- 5. Add 200 μl of Solution P and mix gently by inverting 5 times. Incubate at room temp. for 3 minutes.**

- 6. Centrifuge at 14,000 rpm for 7 minutes. During centrifugation, insert Spin Column into a collection tube.**

- 7. Add 400 μl of Binding Solution to the prepared Spin Column and transfer cleared lysate to the column inserted in the collection tube.**

Avoid disturbing or transferring any of the white precipitate with the supernatant.

- 8. Centrifuge at 10,000 $\times g$ for 30 sec. Discard the flow through and reinsert column in the collection tube.**

Protocol *(continued)*

- 9. Add 700 μ l of Wash Solution . Centrifuge at 10,000 x g for 30 sec. Discard flow through and reinsert column in the collection tube.**
- 10. Centrifuge at 14,000 rpm for 30 sec.**
- 11. Transfer Spin Column to a sterile 1.5 ml microcentrifuge tube. Add 50 μ l of sterile water or a low-ionic strength buffer (10 mM Tris, pH.8.0) to the middle of the Spin Column. Incubate at room temp. for 3 minutes.**

For eluting DNA, use a low-ionic strength buffer (≤ 10 mM in concentration, pH 7-9) or sterile deionized water. For most applications 10 mM Tris base (pH adjusted to 8 with HCl) is recommended; however, TE buffer (10 mM Tris HCl, pH 8.0, 1 mM EDTA) may be used for applications in which EDTA will not interfere with subsequent reactions.

- 12. Centrifuge at 14,000 rpm for 1 minute.**
- 13. Discard the Spin Column, the eluate in the microcentrifuge tube includes plasmid DNA.**

Troubleshooting Guide

| Problem | Reason | Solution |
|---|---|--|
| Poor or low recovery | Improper washing | Confirm the wash solution concentrate was diluted with the specified volume of ethanol. Keep bottles tightly capped between uses to prevent evaporation. |
| | Insufficient number of cells | Culture may be too old. Prepare a new culture. |
| | Prolonged alkaline lysis | Reduce the time for cell lysis to 3 minutes or until the suspended cells form a clear viscous solution after inversion with the lysis solution. |
| | Poor elution | Repeat elution or increase elution volume. Incubation of column at 70°C for 3 minutes with Elution buffer may increase yields. |
| | Lysis reaction may have been incomplete. | Do not exceed recommended culture volumes. Decrease culture volume. |
| Chromosomal DNA contamination | Chromosomal DNA is sheared by vortexing the tube to mix the contents after adding solution L or solution N or solution P. | Mix the contents of the tube by gentle inversion; do not vortex. |
| Enzymatic reactions using recovered DNA do not proceed | DNA concentration is too low. | Precipitate the DNA with alcohol, then resuspend DNA in a smaller volume of elution buffer. |
| | High salt content in the final plasmid DNA eluate | Precipitate the DNA using ethanol. |
| | Residual ethanol from the diluted wash solution | Centrifuge the column for 1 minute after the wash step to remove any residual wash solution. |

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