

HYGİSLİDE TSA / SDA

PRINCIPLE AND INTERPRETATION:

Side1: TSA: A general purpose medium for the growth of a wide variety of organisms. In Tryptic Soy Agar, the combination of casein and soy peptones renders the medium nutritious by supplying organic nitrogen, particularly amino acids and longer-chained peptides. Sodium chloride maintains the osmotic equilibrium. Agar is the solidifying agent.

Side2: SDA: An acidic pH medium for the isolation of dermatophytes, other fungi and yeasts. Sabouraud Dextrose Agar is a peptone medium supplemented with dextrose to support the growth of fungi. The peptones are sources of nitrogenous growth factors. Dextrose provides an energy source for the growth of microorganisms.

COMPOSITION:

TSA

Ingredients	Gr/Liter
Pancreatic digest of casein	15 gr
Enzymatic digest of soya bean	5 gr
Sodium chloride	5 gr
Agar	15 gr

pH: 7,3 ± 0,2

SDA

Ingredients	Gr/Liter
Mycological peptone	10 gr
Glucose(dextrose)	40 gr
Agar	15 gr

pH: 5,6 ± 0,2

***Formula adjusted, standardized to suit performance parameters

INSTRUCTIONS FOR USE:

Testing Fluids:

1. Mix liquid test sample.
2. Remove the paddle from the vial. Do not touch the agar surfaces.
3. Immerse the slide in the fluid to be tested for about 5- 10 seconds. Alternatively expose the slide to a spray or running fluid so that the slide surfaces are covered.
4. Both agar surfaces must be completely contacted.
5. Allow excess fluid to drain off both paddle agar surfaces.
6. Replace the Slide into the tube and twist to tighten the cap. Label the tube with the identification label supplied. Incubate the slide as directed later.

Testing Surfaces:

1. Remove the paddle from the vial. Do not touch the agar surfaces.
2. To assure an accurate area recovery, contact the paddle to 20²cm of the surface by contacting the surface twice in separate 10²cm areas.
3. Replace the Slide into the tube and twist to tighten the cap. Label the tube with the identification label supplied. Incubate the slide as directed later.

QUALITY CONTROL:

1.Sterility Control:

Incubation 2 d at 30-35°C and 3 d at 20-25°C: NO GROWTH

2.Physical/Chemical Control

	pH	Apperance:
TSA:	7,3 ± 0,2	Light amber
SDA:	5,6 ± 0,2	Amber

3.Microbiological Control: Incubate at 35±2 °C temperature for 24 hours and 25±2 °C 48 h-5 days.

Side1: TSA

Microorganism	Inoculum (CFU)	Results	
		Growth	Reaction
Bacillus subtilis ATCC 6633	10-100	Good	Good
Staphylococcus aureus ATCC 6538	10-100	Good	Good
Candida albicans ATCC 10231	10-100	Good	Good
Pseudomonas aeruginosa ATCC 9027	10-100	Good	Good
Aspergillus brasiliensis ATCC 16404	10-100	Good	Good
Escherichia coli ATCC 8739	10-100	Good	Good
Staphylococcus epidermidis ATCC 12228	10-100	Good	Good
Bacillus subtilis ATCC 6633	10-100	Good	Good

Side2: SDA


Microorganism	Inoculum (CFU)	Results	
		Growth	Reaction
Candida albicans ATCC 10231	10-100	Good	Good
Aspegillus brasiliensis ATCC 16404	10-100	Good	Good

INTERPRETATION OF RESULTS

Compare the slide surfaces against the comparison chart printed below. Read the result corresponding to fluids or surfaces as appropriate. Note that very high levels of organisms could lead to a confluent growth and could be recorded as a nil result. Compare against an unused slide when reading results.


Bacteria Comparison Chart

Surfaces
CFU/cm²



Fluids
CFU/mL

Fungi Comparison Chart



DISPOSAL:

Incubated Slides may contain active bacteria and micro-organisms. Do not open infected slides except as part of disposal procedure. Infected slides should be autoclaved, incinerated or opened and soaked in a chlorine-based disinfectant (liquid bleach) for 20 minutes prior to disposal.

STORAGE CONDITIONS AND SHELF LIFE:

Slides should be stored in 2-20 °C. Temperature fluctuations may result in condensation settling at the bottom of the vial, although this does not affect culture properties, it could reduce the shelf-life or cause the agar to separate from the plastic paddle support.

Avoid sudden temperature changes. Shield from direct sunlight. Do not allow paddles to freeze. Do not use any slides which have been inadvertently contaminated during storage and which are already showing growth of micro-organisms

Use before expiry date on the label. Do not use beyond stated expiry date.

PACKAGING:

Katalog Number: 06017

Content/Packaging: 20 Slides/Box

REFERENCES:

1. Abbott J. D. and Graham J. M. (1961) Mon. Bull. Min. Hlth Pub. Hlth Lab. Serv. 20. 51-58.
2. Barrow G. I. and Ellis C. (1962) Mon. Bull. Min. Hlth Pub. Hlth Lab. Serv. 21. 141-147.
3. Cooke G. T. and Daines C. F. (1964) Mon. Bull. Min. Hlth Publ. Hlth Lab. Serv. 23. 81-85.
4. Gillies R. R. (1964) J. Hyg. Camb. 62. 1-9.
5. Mitchell T. G. (1964) J. Appl. Bact. 27. 45-52.
6. Barnes Ella M. and Shrimpton D. H. (1958) J. Appl. Bact. 2. 313-329.
7. American Public Health Association (1978) Standard Methods for the Examination of Dairy Products. 14th Edn. APHA Inc. Washington DC.
8. American Public Health Association, Standard Methods for the Examination of Dairy Products, 14th ed., APHA Inc., Washington, D.C. (1978)
9. E.W. Frampton, et al., Comparison of β -glucuronidase and indole-based direct plating methods for enumeration of unstressed E. coli, J. Food Protect. 53, 933 (1990)
10. Carlier Gwendoline I. M. (1948) Brit. J. Derm. Syph. 60. 61-63.
11. Hodges R. S. (1928) Arch. Derm. Syph., New York, 18. 852.
12. Sabouraud R. (1910) 'Les Teignes', Masson, Paris.
13. Georg Lucille K., Ajello L. and Papageorge Calomira (1954) J. Lab. Clin. Med. 44. 422-428.
14. Ajello Libero (1957) J. Chron. Dis. 5. 545-551.
15. Williams Smith H. and Jones J. E. T. (1963) J. Path. Bact. 86. 387-412.
16. Hantschke D. (1968) Mykosen. 11. 113-115

STERILE A

Aseptic Sterile

LOT

Batch Code

REF

Catalogue Number

CONTROL -

Negative Controls

CONTROL +

Positive Controls



Use by



Temperature
Limitation



Do not reuse



Contains sufficient
for <n> tests



Look at user manual



Manufacturer