



STERIL TRANSPORT SWAB with AMIES AGAR

INTENDED USE:

Transport swab with AMIES is used for collecting and storing microbe specimen. AMIES gel is recommended for most specimens containing aerobic, fastidious and anaerobic bacteria.

COMPOSITION:

Transport swab with AMIES is composed of a swab, a transport tube, tube cap, label and AMIES gel.

AMIES gel Composition:

- Potassium chloride
- Magnesium chloride
- Sodium thioglycollate
- Monobasic potassium phosphate
- Calcium chloride
- Agar

PRECAUTIONS:

- 1 For single use only, discards after use.
- 2 Store in dry and draughty warehouse.
- 3 Please do not use the product if the package is damaged.
4. In vitro diagnostic use only

STERILE CONDITION:

Sterilized by Gamma radiation sterility, the sterilization period is three years.

PACKAGING:

Packaging method: 1pc/blister, 1000pcs/carton.

Katalog Number: 08002

STORAGE CONDITIONS AND SHELF LIFE:

The product should be stored in dry and draughty warehouse. The period of validity is three years since manufacturing date, do not use the product beyond stated expiry date.