



STERIL COTTON SWAB

INTENDED USE:

Cotton swab is mainly used in sampling from the parts of pharyngeal, nasal cavity, auditory meatus and wound.

COMPOSITION:

Cotton swab is composed of wooden stick and tip.

PRECAUTIONS:

- 1 For single use only, discards after use.
- 2 Store in dry and draughty warehouse.
- 3 Please do not use the product if the package is damaged.
4. In vitro diagnostic use only

STERILE CONDITION:

Sterilized by EO gas, the sterilization period is five years.

PACKAGING:

Packaging method: 1pcs/blister, 5000pcs/carton.

Katalog Number: 08005

STORAGE CONDITIONS AND SHELF LIFE:

The product should be stored in dry and draughty warehouse. The period of validity is five years since manufacturing date, do not use the product beyond stated expiry date.