



CERTIFICATE OF ANALYSIS

Product Name	Phosphate Buffered Peptone Saline(PBPS) + %3 Tween80 + %0,3 Lecithin (90mL)
Catalogue No	01152
Batch No	2001115
Dehydrated Brand / Batch No / Expiry Date	Buffered Saline Pepton Water (CONDALAB): 810104 / 10-2022
Neutralizer Brand / Batch No / Expiry Date	Tween 80 (SABO): 20103868 / 05-2020 Lecithin (LAB-M): Q42844 / 10-2020
Production Date	09.01.2020
Expiry Date	09.05.2020
Storage Conditions	2 – 25 °C

Physical & Chemical Specifications

Appearance	clear	Appropriate
Color	Light amber	Appropriate
pH	7,0±0,2	6.86

Microbiological Specifications: Productivity and acc. To ISO/TS 11133

Method of control: Direct Inoculation	Inoculum: :(+ control 10-100 cfu, - control 100-1000 cfu)	Incubation: 35±2 °C:18-24h			
Sterility:	Incubation: (20-25°C : 3 d, 30-35°C : 2 d)	Sterile			
QC Strains	ATCC	Results			
		Recovery		Growth	
<i>Staphylococcus aureus</i>	25923	Good	Good	Good	Appropriate
<i>Pseudomonas aeruginosa</i>	9027	Good	Good	Good	Appropriate
<i>Bacillus subtilis</i>	6633	Good	Good	Good	Appropriate
<i>Escherichia coli</i>	25922	Good	Good	Good	Appropriate
<i>Salmonella typhimurium</i>	14028	Good	Good	Good	Appropriate
<i>Candida albicans</i>	10231	Good	Good	Good	Appropriate

QC Date	14.01.2020
QC Manager	Gizem KOCABAŞ

This document has been prepared electronically.

RTA LABORATUVARLARI BİYOLOJİK ÜRÜNLER İLAÇ VE MAKİNE SANAYİ TİCARET A.Ş.
Gebze Plastikçiler Organize Sanayi Bölgesi (GEPOSB) Cumhuriyet Cad.No.:3 41400 Gebze / KOCAELİ
Tel.:0262 648 53 00 Fax:0262 751 06 77 e-mail: rtalabs.com.tr